

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>07-367-2015</b>	<b>Revision: 5</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Pipette tips PP, premium surface, 0-300µl, 96 pcs/rack, transparent, moulded rings, np pcr ready, autocl. up to 121°C
<b>Dimensions</b>	length: 59,54 mm (+/- 0,2 mm)
	inner diameter: 0,51 mm (+/- 0,2 mm)
<b>Colour</b>	transparent
<b>Graduation</b>	moulded rings
<b>Nominal Capacity</b>	0-300µl
<b>Material Characteristics</b>	polypropylene (PP), premium surface
<b>Autoclavability</b>	autocl. up to 121°C
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile SAL 10-3
<b>Purity Grade</b>	free of RNase and Dnase
	free of pyrogens
	free of PCR inhibitors
	free of human DNA
	free of endotoxins
	free of ATP
<b>Shelf - Life</b>	60 months at controlled storage conditions
<b>Packaging</b>	pcs per inner cs: 960
	pcs per carton: 4.800
	pcs per pallet: 201.600
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

00-000-0001/06.2020



# Technical Datasheet

Date: 03.06.22

Page 2 of 2

Approved by nerbe plus quality department