

Technical Datasheet

| | | |
|----------------------|--------------------|------------------------|
| | | EN ISO 13485 certified |
| Article - No: | 07-369-2015 | Revision: 5 |
| Remarks | | single use only |

Specification

| | |
|---------------------------------|---|
| Product | Pipette tips PP, premium surface, 100-1.250µl, 96 pcs/rack, transparent, moulded rings, np pcr ready, autocl. up to 121°C |
| Dimensions | length: 98,40 mm (+/- 0,2 mm) |
| | inner diameter: 0,89 mm (+/- 0,2 mm) |
| Colour | transparent |
| Graduation | moulded rings |
| Nominal Capacity | 100-1.250µl |
| Material Characteristics | polypropylene (PP), premium surface |
| Autoclavability | autocl. up to 121°C |
| Product Characteristics | biologically inert |
| | physiologically harmless |
| | no heavy metals according to EC/94/62 |
| | sterile SAL 10-3 |
| Purity Grade | free of RNase and Dnase |
| | free of pyrogens |
| | free of PCR inhibitors |
| | free of human DNA |
| | free of endotoxins |
| | free of ATP |
| Shelf - Life | 60 months at controlled storage conditions |
| Packaging | pcs per inner cs: 768 |
| | pcs per carton: 3.840 |
| | pcs per pallet: 107.520 |
| Storage | dry condition: 40%-60% rel. humidity |
| | bags / cases orig. closed |
| | storage temp. +10°C to max. +30°C |
| | pal. not stackable |

Technical Datasheet

Date: 03.06.22

Page 2 of 2

Approved by nerbe plus quality department