

Technical Datasheet

		EN ISO 13485 certified
Article - No:	07-375-2015	Revision: 6
Remarks		single use only

Specification

Product	Pipette tips PP, premium surface, refill system, 0-200µl, 96 pcs/tray, transparent tips, translucent yellow tray, grad., np pcr ready, autocl. up to 121°C
Dimensions	length: 50,80 mm (+/- 0,2 mm)
	inner diameter: 0,53 mm (+/- 0,2 mm)
Colour	transparent tips, translucent yellow tray
Graduation	grad.
Nominal Capacity	0-200µl
Material Characteristics	polypropylene (PP), premium surface
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile SAL 10-3
Purity Grade	free of RNase and Dnase
	free of pyrogens
	free of PCR inhibitors
	free of human DNA
	free of endotoxins
	free of ATP
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 960
	pcs per carton: 4.800
	pcs per pallet: 470.400
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable



Technical Datasheet

Date: 03.06.22

Page 2 of 2

Approved by nerbe plus quality department