

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>07-622-0000</b>	<b>Revision: 5</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Filter tip PP, premium surface, 0-20µl, transparent, moulded rings, np pcr ready
<b>Dimensions</b>	length: 50,80 mm (+/- 0,2 mm)
	inner diameter: 0,53 mm (+/- 0,2 mm)
<b>Colour</b>	transparent
<b>Graduation</b>	moulded rings
<b>Nominal Capacity</b>	0-20µl
<b>Material Characteristics</b>	polypropylene (PP), premium surface
	filter: high density polyethylene (PE)
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile SAL 10-3
<b>Purity Grade</b>	free of RNase and Dnase
	free of pyrogens
	free of PCR inhibitors
	free of human DNA
	free of endotoxins
	free of ATP
<b>Shelf - Life</b>	60 months at controlled storage conditions
<b>Packaging</b>	pcs per bag: 1.000
	pcs per carton: 10.000
	pcs per pallet: 280.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

00-000-0001/06.2020



# Technical Datasheet

Date: 03.06.22

Page 2 of 2



EN ISO 13485 certified

Approved by nerbe plus quality department

00-000-0001/06.2020

nerbe plus GmbH & Co. KG  
Porschestr. 25  
D-21423 Winsen/Luhe  
Tel: +49 (0)4171 7849 0  
Fax: +49 (0)4171 6 45 45

info@nerbe-plus.de  
www.nerbe-plus.de  
HRA Lüneburg 203240  
USt.-ID: DE 32583 0652

Pers. haftende Gesellschafterin:  
NERBE R&D GmbH, Winsen/Luhe  
GF: Philipp J.F. Nerbe  
HRB Lüneburg 208826

Deutsche Bank AG  
Kto. 0204 41200  
IBAN: DE41 2407 0075 0020 4412 00  
SWIFT / BIC DEUTDE2H240

Hamburger Sparkasse  
Kto. 1387 140831  
IBAN: DE56 2005 0550 1387 1408 31  
SWIFT / BIC HASPDE33XXX

