

Technical Datasheet

| | | |
|----------------------|--------------------|------------------------|
| | | EN ISO 13485 certified |
| Article - No: | 07-642-0000 | Revision: 5 |
| Remarks | | single use only |

Specification

| | |
|---------------------------------|---|
| Product | Filter tip PP, premium surface, 0-100µl, transparent, moulded rings, np pcr ready |
| Dimensions | length: 59,54 mm (+/- 0,2 mm) |
| | inner diameter: 0,51 mm (+/- 0,2 mm) |
| Colour | transparent |
| Graduation | moulded rings |
| Nominal Capacity | 0-100µl |
| Material Characteristics | polypropylene (PP), premium surface |
| | filter: high density polyethylene (PE) |
| Product Characteristics | biologically inert |
| | physiologically harmless |
| | no heavy metals according to EC/94/62 |
| | sterile SAL 10-3 |
| Purity Grade | free of RNase and Dnase |
| | free of pyrogens |
| | free of PCR inhibitors |
| | free of human DNA |
| | free of endotoxins |
| | free of ATP |
| Shelf - Life | 60 months at controlled storage conditions |
| Packaging | pcs per bag: 1.000 |
| | pcs per carton: 10.000 |
| | pcs per pallet: 280.000 |
| Storage | dry condition: 40%-60% rel. humidity |
| | bags / cases orig. closed |
| | storage temp. +10°C to max. +30°C |
| | pal. not stackable |

Technical Datasheet

Date: 03.06.22

Page 2 of 2



EN ISO 13485 certified

Approved by nerbe plus quality department

nerbe plus GmbH & Co. KG
Porschestr. 25
D-21423 Winsen/Luhe
Tel: +49 (0)4171 7849 0
Fax: +49 (0)4171 6 45 45

info@nerbe-plus.de
www.nerbe-plus.de
HRA Lüneburg 203240
USt.-ID: DE 32583 0652

Pers. haftende Gesellschafterin:
NERBE R&D GmbH, Winsen/Luhe
GF: Philipp J.F. Nerbe
HRB Lüneburg 208826

Deutsche Bank AG
Kto. 0204 41200
IBAN: DE41 2407 0075 0020 4412 00
SWIFT / BIC DEUTDE2H240

Hamburger Sparkasse
Kto. 1387 140831
IBAN: DE56 2005 0550 1387 1408 31
SWIFT / BIC HASPDE33XXX

