

Technical Datasheet

EN ISO 13485 certified		
Article - No:	07-895-0500	Revision: 5
Remarks	suitable for: Gilson P 5000 and F 123603, MLA (universal), Eppendorf (universal)	single use only

Specification

Product	Pipette tip PP "classic", 10.000µl, universal fit (Gilson, Eppendorf, Labnet), transparent, grad., np pcr ready, autocl. up to 121°C
Dimensions	length: 243,00 mm (+/- 0,2 mm)
Colour	transparent
Graduation	grad.
Nominal Capacity	10.000µl
Material Characteristics	polypropylene (PP)
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 100
	pcs per carton: 1.000
	pcs per pallet: 20.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department

