

Technical Datasheet

		EN ISO 13485 certified
Article - No:	08-002-6004	Revision: 2
Remarks		single use only

Specification

Product	Universal specimen container with spoon PP, 12ml, Ø16x100 mm, with fitted screw cap PE, container transparent, cap brown, grad., sterile R, max. RCF 2.500g
Dimensions	outer diameter: 16,00 mm (+/- 0,2 mm)
	height (without screw cap): 100,00 mm (+/- 0,2 mm)
Colour	container transparent, cap brown
Graduation	grad.
Nominal Capacity	12ml
Material Characteristics	tube: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-20°C - 65°C
Centrifugation	max. RCF 2.500g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per bag: 100
	pcs per carton: 1.000
	pcs per pallet: 24.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department

