

# Technical Datasheet

Date: 22.06.20

Page 1 of 1

EN ISO 13485 certified		
Article - No:	08-122-2004	Revision: 1
Remarks		single use only

## Specification

Product	Universal specimen container with spoon PP, 30ml, Ø25x107 mm, with fitted screw cap PE, skirted, container transparent, cap brown, graduated, max. RCF 2.500g, CE/IVD
Classification	CE/IVD
Dimensions	diameter: 25,00 mm (+/- 0,2 mm)
	height (without screw cap): 107,00 mm (+/- 0,2 mm)
Colour	container transparent, cap brown
Graduation	grad.
Nominal Capacity	30ml
Material Characteristics	tube: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-20°C - 65°C
Centrifugation	max. RCF 2.500g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 100
	pcs per carton: 500
	pcs per pallet: 12.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	