

Technical Datasheet

		EN ISO 13485 certified
Article - No:	08-222-1000	Revision: 2
Remarks		single use only

Specification

Product	Universal specimen container PP, 15ml, Ø20x76mm, with enclosed screw cap PE, skirted, tube transparent, cap natural, grad., max. RCF 4.000g
Dimensions	diameter: 20,00 mm (+/- 0,2 mm)
	height (without screw cap): 76,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap natural
Graduation	grad.
Nominal Capacity	15ml
Material Characteristics	container: polypropylene (PP)
	cap: polyethylene (PE)
Centrifugation	max. RCF 4.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 1.000
	pcs per carton: 1.000
	pcs per pallet: 24.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	