

Technical Datasheet

		EN ISO 13485 certified
Article - No:	08-400-6000	Revision: 2
Remarks		single use only

Specification

Product	Sputum complete transport system, 24ml, content: inner container PS/PE 08-411-6000 and outer container PP/PE 08-422-1000
Dimensions	diameter: 79,00 mm (+/- 0,2 mm)
	height: 28,00 mm (+/- 0,2 mm)
Nominal Capacity	24ml
Material Characteristics	container: polysterene (PS)
	plug: polyethylene (PE)
Temperature range	-20°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Sterilisation Process	radiation sterilization SAL 10-6 inner container and lid only
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 100
	pcs per carton: 100
	pcs per pallet: 4.800
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

