

Technical Datasheet

		EN ISO 13485 certified
Article - No:	08-422-1000	Revision: 2
Remarks		single use only

Specification

Product	Outer container PP for sputum complete system, 24ml, Ø79x28mm, with enclosed lid PE, container natural, lid yellow
Dimensions	diameter: 79,00 mm (+/- 0,2 mm)
	height: 28,00 mm (+/- 0,2 mm)
Colour	container natural, lid yellow
Nominal Capacity	24ml
Material Characteristics	container: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-10°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 100
	pcs per carton: 100
	pcs per pallet: 4.800
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

