

Technical Datasheet

		EN ISO 13485 certified
Article - No:	08-602-0030	Revision: 4
Remarks		single use only

Specification

Product	Universal container PP, 125ml, Ø55x70 mm, with screw thread, without cap, transparent, grad., writ. area, np pcr ready, autocl. up to 121°C
Dimensions	upper diameter: 48,00 mm (+/- 0,2 mm)
	bottom diameter: 55,00 mm (+/- 0,2 mm)
Colour	transparent
Graduation	grad., writ. area
Nominal Capacity	125ml
Material Characteristics	polypropylene (PP)
Temperature range	-40°C-75°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Packaging	pcs per bag: 750
	pcs per carton: 750
	pcs per pallet: 18.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	