

Technical Datasheet

		EN ISO 13485 certified
Article - No:	08-602-0040	Revision: 1
Remarks		single use only

Specification

Product	Screw cap PE, suitable for universal container 125ml, red, np pcr ready
Dimensions	diameter: 55,00 mm (+/- 0,2 mm)
	height: 15,00 mm (+/- 0,2 mm)
Colour	red
Material Characteristics	polyethylene (PE)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Packaging	pcs per bag: 500
	pcs per carton: 1.000
	pcs per pallet: 24.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	