

Technical Datasheet

EN ISO 13485 certified	
Article - No:	08-602-1903
Revision:	2
Remarks	single use only

Specification

Product	Universal container PP, 125ml, Ø55x72 mm, with enclosed screw cap PE, container transparent, cap red, graduated, with writing area, np pcr ready, Autoclaveable up to 121°C (no cap), CE/IVD
Classification	CE/IVD
Dimensions	bottom Ø: 48,0mm (+/- 0,2mm)
	upper Ø: 55,0mm (+/- 0,2mm)
	height (w/o cap): 71,7mm (+/- 0,2mm)
	height (w/cap): 72,0mm (+/- 0,2mm)
Colour	container transp., cap red
Graduation	grad., writ. area
Nominal Capacity	125ml
Material Characteristics	cont.: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-40°C-75°C
Autoclavability	Autocl. up to 121°C (no cap)
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Packaging	pcs per bag: 500
	pcs per carton: 500
	pcs per pallet: 12.000
Storage	dry condition: 40%-60% rel. humidity



Technical Datasheet

Date: 22.06.20

Page 2 of 2

	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	