

Technical Datasheet

		EN ISO 13485 certified
Article - No:	08-602-6041	Revision: 4
Remarks		single use only

Specification

Product	Universal container PP, 125ml, Ø55x72 mm, with fitted screw cap PE, container transp., cap blue grad., writ. area, np pcr ready, sterile , Autocl. up to 121°C (no cap)
Dimensions	bottom Ø: 48,0mm (+/- 0,2mm)
	upper Ø: 55,0mm (+/- 0,2mm)
	height (w/o cap): 71,7mm (+/- 0,2mm)
	height (w/cap): 72,0mm (+/- 0,2mm)
Colour	container transp., cap blue
Graduation	grad., writ. area
Nominal Capacity	125ml
Material Characteristics	cont.: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-40°C-75°C
Autoclavability	Autocl. up to 121°C (no cap)
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per bag: 250
	pcs per carton: 250
	pcs per pallet: 5.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed

00-000-0001/06.2020



Technical Datasheet

	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

