

Technical Datasheet

		EN ISO 13485 certified
Article - No:	08-602-7803	Revision: 3
Remarks		single use only

Specification

Product	Universal container PP, 125ml, Ø55x72mm, fitted screw cap PE, label, seal of originality, container transp., cap yellow, grad., writ. area, np pcr ready, sterile
Dimensions	bottom Ø: 48,0mm (+/- 0,2mm)
	upper Ø: 55,0mm (+/- 0,2mm)
	height (w/o cap): 71,7mm (+/- 0,2mm)
	height (w/cap): 72,0mm (+/- 0,2mm)
Colour	container transp., cap yellow
Graduation	grad., writ. area
Nominal Capacity	125ml
Material Characteristics	cont.: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-40°C-75°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per bag: 240
	pcs per carton: 240
	pcs per pallet: 5.760
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C



Technical Datasheet

Date: 10.06.22

Page 2 of 2

pal. not stackable

Approved by nerbe plus quality department
--