

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>08-622-0000</b>	<b>Revision: 1</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Universal container PP, 125ml, Ø50x68 mm, without lid, transparent, graduated, with writing area, autoclavable up to 121°C, CE/IVD
<b>Classification</b>	CE/IVD
<b>Dimensions</b>	diameter on the bottom: 50,00 mm (+/- 0,2mm) height: 68,00 mm (+/- 0,2mm)
<b>Colour</b>	transparent
<b>Graduation</b>	grad., writ. area
<b>Nominal Capacity</b>	125ml
<b>Material Characteristics</b>	polypropylene (PP)
<b>Temperature range</b>	-196°C - 121°C
<b>Autoclavability</b>	autocl. up to 121°C
<b>Product Characteristics</b>	biologically inert physiologically harmless no heavy metals according to EC/94/62
<b>Packaging</b>	pcs per bag: 50 pcs per carton: 1.000 pcs per pallet: 24.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
<b>Approved by nerbe plus quality department</b>	

00-000-0001/06.2020

