

Technical Datasheet

Date: 22.06.20

Page 1 of 1

		EN ISO 13485 certified
Article - No:	08-622-1902	Revision: 1
Remarks		single use only

Specification

Product	Universal container PP, 125ml, Ø50x68 mm, with lid PE with spout, container transparent, lid red graduated, with writing area, CE/IVD
Classification	CE/IVD
Dimensions	diameter on the bottom: 50,00 mm (+/- 0,2mm) height: 68,00 mm (+/- 0,2mm)
Colour	container transparent, lid red
Graduation	grad., writ. area
Nominal Capacity	125ml
Material Characteristics	container: polypropylene (PP) cap: polyethylene (PE)
Temperature range	-50°C - 60°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Packaging	pcs per bag: 50 pcs per carton: 500 pcs per pallet: 12.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	