

# Technical Datasheet

Date: 22.06.20

Page 1 of 1

		EN ISO 13485 certified
<b>Article - No:</b>	<b>08-633-0122</b>	<b>Revision: 1</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Lid PE, with spout, suitable for universal container 08-622-0000, red
<b>Colour</b>	red
<b>Material Characteristics</b>	polyethylene (PE)
<b>Temperature range</b>	-50°C - 60°C
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
<b>Packaging</b>	pcs per bag: 500
	pcs per carton: 1.000
	pcs per pallet: 24.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	

