

Technical Datasheet

| | |
|------------------------|------------------------|
| EN ISO 13485 certified | |
| Article - No: | 08-733-3004 |
| Revision: | 1 |
| Remarks | single use only |

Specification

| | |
|---------------------------------|--|
| Product | 24-hrs-urine collection bottle PP, with fitted screw cap PE, 2,7 liter, with insert, bottle brown, screw cap yellow, graduated, CE/IVD |
| Classification | CE/IVD |
| Dimensions | length: 112,00 mm (+/- 0,2 mm) |
| | width: 117,00 mm (+/- 0,2 mm) |
| | height: 250,00 mm (+/- 0,2 mm) |
| Colour | bottle brown, screw cap yellow |
| Graduation | grad. |
| Nominal Capacity | 2,7 liter |
| Material Characteristics | bottle: polypropylene (PP) |
| | screw cap: polyethylene (PE) |
| Temperature range | -40°C - 120°C |
| Product Characteristics | biologically inert |
| | physiologically harmless |
| | no heavy metals according to EC/94/62 |
| | Food contact: raw material in compliance with regulation EC 1935/2004 |
| Packaging | pcs per cs: 30 |
| | pcs per carton: 30 |
| | pcs per pallet: 360 |
| Storage | dry condition: 40%-60% rel. humidity |
| | bags / cases orig. closed |
| | storage temp. +10°C to max. +30°C |
| | pal. not stackable |

Approved by nerbe plus quality department

