

Technical Datasheet

		EN ISO 13485 certified
Article - No:	08-733-3005	Revision: 1
Remarks		single use only

Specification

Product	24-hrs-urine collection bottle PP, with fitted screw cap PE, 2,7 liter, with label, with insert bottle brown, screw cap yellow, graduated, CE/IVD
Classification	CE/IVD
Dimensions	length: 112,00 mm (+/- 0,2 mm) width: 117,00 mm (+/- 0,2 mm) height: 250,00 mm (+/- 0,2 mm)
Colour	bottle brown, screw cap yellow
Graduation	grad.
Nominal Capacity	2,7 liter
Material Characteristics	bottle: polypropylene (PP) screw cap: polyethylene (PE)
Temperature range	-40°C - 120°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 Food contact: raw material in compliance with regulation EC 1935/2004
Packaging	pcs per carton: 30 pcs per pallet: 360
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

