

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-001-0010	Revision: 3
Remarks		single use only

Specification

Product	Petri dish PS, Ø55x14,2 mm, without vents, transparent
Dimensions	diameter: 55,00 mm (+/- 0,2 mm)
	height: 14,20 mm (+/- 0,2 mm)
Colour	transparent
Material Characteristics	polystyrene (PS)
Temperature range	-20°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	sterile SAL 10-3
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 15
	pcs per carton: 1.005
	pcs per pallet: 32.160
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	