

Technical Datasheet

EN ISO 13485 certified		
Article - No:	09-021-0010	Revision: 4
Remarks	max. volume 80ml, working volume 15 ml	single use only

Specification

Product	Petri dish PS, Ø92x16 mm, without vents, (non slippery / stackable), transparent
Dimensions	diameter: 92,00 mm (+/- 0,2 mm)
	height: 16,00 mm (+/- 0,2 mm)
Colour	transparent
Material Characteristics	polystyrene (PS)
Temperature range	-20°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile SAL 10-3
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 20
	pcs per carton: 480
	pcs per pallet: 9.600
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

