

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-031-0040	Revision: 4
Remarks		single use only

Specification

Product	Petri dish PS, Ø91,5x14,55 mm, with 3 vents, (non slippery / stackable), transparent
Dimensions	diameter: 91,50 mm (+/- 0,2 mm)
	height: 14,55 mm (+/- 0,2 mm)
Colour	transparent
Material Characteristics	polystyrene (PS)
Temperature range	-20°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 20
	pcs per carton: 480
	pcs per pallet: 9.600
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	