

Technical Datasheet

Date: 22.06.20

Page 1 of 1

		EN ISO 13485 certified
Article - No:	09-061-5010	Revision: 3
Remarks		single use only

Specification

Product	Petri dish PS, Ø140x20 mm, with 3 vents, transparent, sterile R
Dimensions	diameter: 140,00 mm (+/- 0,2 mm)
	height: 20,00 mm (+/- 0,2 mm)
Colour	transparent
Material Characteristics	polystyrene (PS)
Temperature range	-20°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile R SAL 10-6
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 10
	pcs per carton: 110
	pcs per pallet: 3.520
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	