

Technical Datasheet

EN ISO 13485 certified		
Article - No:	09-081-0010	Revision: 4
Remarks	Moulded grid 10 x 10	single use only

Specification

Product	Contact plate PS, Ø65x14,5 mm, with 4 vents, transparent
Dimensions	diameter: 65,00 mm (+/- 0,2 mm)
	height: 14,50 mm (+/- 0,2 mm)
Colour	transparent
Material Characteristics	polystyrene (PS)
Temperature range	-20°C - 50°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 18
	pcs per carton: 918
	pcs per pallet: 29.376
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	