

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-101-0051	Revision: 3
Remarks		single use only

Specification

Product	Petri dish PS, Ø92x14,2 mm, with 3 vents, half division, transparent, sterile
Dimensions	diameter: 92,00 mm (+/- 0,2 mm)
	height: 14,20 mm (+/- 0,2 mm)
Colour	transparent
Material Characteristics	polystyrene (PS)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 25
	pcs per carton: 600
	pcs per pallet: 12.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	