

# Technical Datasheet

EN ISO 13485 certified		
<b>Article - No:</b>	<b>09-302-0020</b>	<b>Revision: 2</b>
<b>Remarks</b>	<b>max. vol.: approx. 3 l / working vol.: approx. 1,25 l</b>	<b>single use only</b>

## Specification

<b>Product</b>	Autoclavable bag CPP, 200x300 mm, thickness: 50µm natural, autoclavable up to 134°C
<b>Dimensions</b>	width: 200,00 mm (+/- 6 mm)
	length: 300,00 mm (+/- 15 mm)
	thickness: 50 µm (+/- 10%)
<b>Colour</b>	natural
<b>Material Characteristics</b>	cast polypropylene (CPP)
<b>Autoclavability</b>	autocl. up to 134°C
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	side seam
	density approx. 0,90 gr.
<b>Packaging</b>	pcs per carton: 1.000
	pcs per pallet: 60.000
	pcs per dsp case: 100
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	

