

Technical Datasheet

Date: 22.06.20

Page 1 of 1

EN ISO 13485 certified		
Article - No:	09-312-0000	Revision: 2
Remarks	max. vol.: approx. 11 l / working vol.: approx. 5 l	single use only

Specification

Product	Autoclavable bag PP, 300x500 mm, thickness: 50µm, natural, autoclavable up to 134°C
Dimensions	width: 300,00 mm (+/- 8 mm)
	length: 500,00 mm (-10/+20 mm)
	thickness: 50 µm (+/- 10%)
Colour	natural
Material Characteristics	polypropylene (PP)
Autoclavability	autocl. up to 134°C
Product Characteristics	biologically inert
	physiologically harmless
	extremely resilient double bottom seam
	density approx. 0,90 gr.
Packaging	pcs per cs: 1.000
	pcs per carton: 1.000
	pcs per pallet: 20.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	