

Technical Datasheet

Date: 22.06.20

Page 1 of 1

EN ISO 13485 certified		
Article - No:	09-372-0070	Revision: 2
Remarks	max. vol.: approx. 120 l / working vol.: approx. 60 l	single use only

Specification

Product	Autoclavable bag PP, 700x1.100 mm, thickness: 70µm natural, autoclavable up to 134°C
Dimensions	width: 700,00 mm (+/- 8 mm)
	length: 1.100,00 mm (-10/+20 mm)
	thickness: 70 µm (+/- 10%)
Colour	natural
Material Characteristics	polypropylene (PP)
Autoclavability	autocl. up to 134°C
Product Characteristics	biologically inert
	physiologically harmless
	extremely resilient double bottom seam
	density approx. 0,90 gr.
Packaging	pcs per cs: 300
	pcs per carton: 300
	pcs per pallet: 5.400
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	