

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>09-401-5003</b>	<b>Revision: 4</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Inoculating loop PS, 1µl, flexible, 20 pcs Zip Bag dark green, sterile R
<b>Dimensions</b>	diameter loop: 0,90 mm (+/- 0,1 mm)
	diameter handle: 4,00 mm (+/- 0,1 mm)
	length: 199,00 mm (+/- 0,2 mm)
<b>Colour</b>	dark green
<b>Nominal Capacity</b>	1µl
<b>Material Characteristics</b>	polystyrene (PS)
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
<b>Purity Grade</b>	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
<b>Sterilisation Process</b>	radiation sterilization SAL 10-6
<b>Shelf - Life</b>	36 months at controlled storage conditions
<b>Packaging</b>	pcs per inner cs: 1.000
	pcs per carton: 5.000
	pcs per pallet: 180.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	

