

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-420-5063	Revision: 6
Remarks		single use only

Specification

Product	Inoculating loop PS, 10µl, rigid, 20 pcs Zip Bag, transparent, color code: turquoise, sterile R
Colour	transparent, color code: turquoise
Nominal Capacity	10µl
Material Characteristics	polystyrene (PS)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 1.000
	pcs per carton: 5.000
	pcs per pallet: 180.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

