

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-421-5002	Revision: 2
Remarks		single use only

Specification

Product	Inoculating loop/needle PS, 10µl, flexible, Ø4x200 mm, 20pcs/peel-pack, light blue, np pcr ready, sterile EO
Dimensions	diameter loop: 3,30 mm (+/- 0,1 mm)
	diameter handle: 3,00 mm (+/- 0,1 mm)
	length: 192,00 mm (+/- 0,2 mm)
Colour	light blue
Nominal Capacity	10µl
Material Characteristics	polystyrene (PS)
Temperature range	-50°C - 55°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	gas sterilisation by ethyleneoxid
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 1.000
	pcs per carton: 5.000
	pcs per pallet: 180.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	