

# Technical Datasheet

Date: 22.06.20

Page 1 of 1

		EN ISO 13485 certified
<b>Article - No:</b>	<b>09-453-0100</b>	<b>Revision: 1</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Seal bags LDPE, 200x3000 mm, thickness: 50µm, no printing, with PE-Logo, transparent
<b>Dimensions</b>	width: 200,00 mm (+/- 6 mm)
	length: 300,00 mm (+/- 15 mm)
	thickness: 50 µm (+/- 10%)
<b>Colour</b>	transparent
<b>Material Characteristics</b>	low density polyethylene (LDPE)
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
<b>Packaging</b>	pcs per bag: 100
	pcs per carton: 1.000
	pcs per pallet: 64.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	