

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-511-5017	Revision: 8
Remarks		single use only

Specification

Product	Culture swab, without transport medium, flexible PS stick with viscose tip, color code white, sterile R, CE/IVD
Classification	CE/IVD
Dimensions	tube:
	length: 152,00 mm (+/- 5 mm)
	outer diameter: 11,50 mm (+/- 0,5 mm)
	inner diameter: 10,00 mm (+/- 0,5 mm)
	plug:
	length: 26,00 mm (+/- 0,2 mm)
	outer diameter: 12,00 mm (+/- 0,5 mm)
	stick:
	length: 150,00 mm (+/- 3 mm)
	diameter: 2,50 mm (+/- 0,2 mm)
diameter tip: 5,00 mm (+/- 1 mm)	
total length of tube and plug: 174,00 mm (+/- 3 mm)	
Material Characteristics	tube: polypropylene (PP)
	plug: polyethylene (PE)
	swab: polysterene (PS)
	tip of swab: viscose (VI)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	free of detectable phthalate acc. to REACH regulation
	free of latex
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Sterilisation Process	radiation sterilization SAL 10-6



Technical Datasheet

Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per carton: 2.000
	pcs per pallet: 48.000
	pcs per dsp case: 500
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	