

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-551-8066	Revision: 3
Remarks		single use only

Specification

Product	Transport swab with amies medium, flexible PS stick with viscose tip, color code blue, single peel-pack, sterile R, CE/IVD
Classification	CE/IVD
Dimensions	tube:
	length: 145,00 mm (+/- 0,2 mm)
	outer diameter: 12,00 mm (+/- 0,2 mm)
	inner diameter: 10,70 mm (+/- 0,2 mm)
	plug:
	length: 20,00 mm (+/- 0,2 mm)
	outer diameter: 12,00 mm (+/- 0,2 mm)
	stick:
	length: 150,00 mm (+/- 0,2 mm)
	diameter: 2,50 mm (+/- 0,2 mm)
	diameter tip: 4,00 mm (+/- 0,2 mm)
	peelpack:
length: 38,00 mm (+/- 0,2 mm)	
width: 210,00 mm (+/- 0,2 mm)	
total length of tube and plug: 165,00 mm (+/- 0,2 mm)	
Material Characteristics	tube: polypropylene (PP)
	plug: polyethylene (PE)
	swab: polysterene (PS)
	tip of swab: viscose (VI)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	free of detectable phthalate acc. to REACH regulation



Technical Datasheet

	free of latex
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	30 months at controlled storage conditions
Packaging	pcs per carton: 1.200
	pcs per pallet: 19.200
	pcs per dsp case: 150
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	