

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-580-0110	Revision: 6
Remarks		single use only

Specification

Product	, 1 flocked swab (oral) with nylon head, sterile R, breakpoint at 30mm and 80mm and ABS stick, , tube CE/IVD, swab CE2797
Classification	tube CE/IVD, swab CE2797
Dimensions	Tube:
	Length: 82,6 mm
	Diameter: 13 mm
	Cap Diameter: 16,6 mm
	Swab (rod):
	Length: 151 mm (+/-0,3mm)
	Diameter: 2,55 mm und 1,75 mm (+/-0,15mm)
	Tip Diameter: 3,5 mm (+/-1,15mm)
Tip Length: 16,5 mm (+4mm)	
Breakpoint: 80 mm	
Material Characteristics	tube: polypropylene (PP)
	tube cap: polyethylene (PE)
	swab: Acrylonitrile butadiene styrene (ABS)
	head of swab: nylon
Temperature range	+5 °C - +25°C
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	medium: 18 months / swab: 36 months at controlled storage conditions
Packaging	pcs per inner cs: 100
	pcs per carton: 600
	pcs per pallet: 19.200
Storage	dry condition: 40%-60% rel. humidity
	inner box / cases orig. closed
	storage temp. +5°C to max. +25°C

Technical Datasheet

	pal. not stackable
Approved by nerbe plus quality department	

