

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-580-0120	Revision: 6
Remarks		single use only

Specification

Product	np GREEN SWAB, LTM, 3ml , 1 flocked swab (nasall) with nylon head, sterile R, breakpoint at 80mm and ABS stick, , , , tube CE/IVD, swab CE2797
Classification	tube CE/IVD, swab CE2797
Dimensions	Tube:
	Length: 82,6 mm
	Diameter: 13 mm
	Cap Diameter: 16,6 mm
	Swab (rod):
	Length: 150 mm
	Diameter: 2,5 mm and 1,5 mm
	Tip Diameter: 2mm
Tip Length: 20mm	
Breakpoint: 80 mm	
Material Characteristics	tube: polypropylene (PP)
	tube cap: polyethylene (PE)
	swab: Acrylonitrile butadiene styrene (ABS)
	head of swab: nylon
Temperature range	+5 °C - +25°C
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	medium: 18 months / swab: 36 months at controlled storage conditions
Packaging	pcs per inner cs: 100
	pcs per carton: 600
	pcs per pallet: 19.200
Storage	dry condition: 40%-60% rel. humidity
	inner box / cases orig. closed
	storage temp. +5°C to max. +25°C



Technical Datasheet

Date: 24.09.21

Page 2 of 2

pal. not stackable

Approved by nerbe plus quality department
--