

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-613-5001	Revision: 6
Remarks		single use only

Specification

Product	Blender bags PE, 400ml, 180x300 mm, without filter, 25 pcs/bag, transparent, np pcr ready, sterile
Dimensions	width: 180,00 mm (+/- 0,2 mm)
	length: 300,00 mm (+/- 0,2 mm)
	thickness: 70 µm
Colour	transparent
Nominal Capacity	400ml
Material Characteristics	polyethylene (PE)
Temperature range	-40°C - 80°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	10 years under controlled storage condition
Packaging	pcs per inner cs: 500
	pcs per carton: 2.500
	pcs per pallet: 60.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

