

Technical Datasheet

EN ISO 13485 certified		
Article - No:	09-681-0000	Revision: 4
Remarks	Suitable bags: 09-613-5001, 09-613-5031, 09-613-5041	single use only

Specification

Product	Closing-Clip PVC, for blender bags 400ml, 50 pcs/bag, blue
Colour	blue
Material Characteristics	Polyvinyl chloride (PVC)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Packaging	pcs per inner cs: 200
	pcs per carton: 1.000
	pcs per pallet: 32.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	