

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-723-5000	Revision: 1
Remarks		single use only

Specification

Product	Wire closure bags PE, 1.260ml, 140x380 mm, transparent, sterile R
Dimensions	width: 140,00 mm (+/- 0,2 mm)
	length: 380,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	1.260ml
Material Characteristics	bag: polyethylene (PE)
	closure: steel wires covered with pressure sensitive tape with a backing of
	polypropylene (PP) or Polyvinyl chloride (PVC)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	leak-proof
	odor-proof
Purity Grade	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 500
	pcs per carton: 1.000
	pcs per pallet: 65.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

