

Technical Datasheet

EN ISO 13485 certified		
Article - No:	09-733-5000	Revision: 1
Remarks	material thickness: 0,08 mm	single use only

Specification

Product	Wire closure bags PE, 1.650ml, 178x305 mm, transparent, sterile R
Dimensions	width: 178,00 mm (+/- 0,2 mm)
	length: 305,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	1.650ml
Material Characteristics	bag: polyethylene (PE)
	closure: steel wires covered with pressure sensitive tape with a backing of polypropylene (PP) or Polyvinyl chloride (PVC)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	leak-proof
	odor-proof
Purity Grade	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 250
	pcs per carton: 1.000
	pcs per pallet: 52.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

