

Technical Datasheet

| | | |
|----------------------|--------------------|------------------------|
| | | EN ISO 13485 certified |
| Article - No: | 09-812-0001 | Revision: 5 |
| Remarks | | single use only |

Specification

| | |
|--|---|
| Product | PS stick with small viscose tip, Ø2,5x150mm, CE |
| Classification | CE |
| Dimensions | diameter: 2,50 mm (+/- 0,2 mm) |
| | length: 150,00 mm (+/- 0,2 mm) |
| | diameter tip : 4,00-5,50 mm (+/- 0,2 mm) |
| Material Characteristics | stick: polystyrene (PS) |
| | tip: viscose (VI) |
| Product Characteristics | biologically inert |
| | physiologically harmless |
| | free of detectable phthalate acc. to REACH regulation |
| Purity Grade | manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1) |
| Sterilisation Process | N/A |
| Shelf - Life | 60 months at controlled storage conditions |
| Packaging | pcs per bag: 100 |
| | pcs per carton: 9.000 |
| | pcs per pallet: 504.000 |
| Storage | dry condition: 40%-60% rel. humidity |
| | bags / cases orig. closed |
| | storage temp. +10°C to max. +30°C |
| | pal. not stackable |
| Approved by nerbe plus quality department | |

