

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-812-8052	Revision: 7
Remarks		single use only

Specification

Product	PS stick with small viscose tip, Ø2,5x150mm, single peel-pack, sterile EO, CE2797
Classification	CE2797
Dimensions	diameter: 2,50 mm (+/- 0,2 mm)
	length: 150,00 mm (+/- 0,2 mm)
	diameter tip : 4,00-5,50 mm (+/- 0,2 mm)
Material Characteristics	stick: polystyrene (PS)
	tip: viscose (VI)
Product Characteristics	biologically inert
	physiologically harmless
	free of detectable phthalate acc. to REACH regulation
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Sterilisation Process	gas sterilisation by ethyleneoxid, 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 100
	pcs per carton: 1.000
	pcs per pallet: 70.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

