

# Technical Datasheet

EN ISO 13485 certified		
<b>Article - No:</b>	<b>09-812-8054</b>	<b>Revision: 8</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	PS stick with breakpoint, small Viscose tip, Ø2,5x150mm, single peel-pack, sterile EO, CE2797
<b>Classification</b>	CE2797
<b>Dimensions</b>	diameter: 2,50 mm (+/- 0,2 mm) length: 150,00 mm (+/- 0,2 mm) diameter tip : 4,00-5,50 mm (+/- 0,2 mm)
<b>Material Characteristics</b>	stick: polystyrene (PS) tip: viscose (VI)
<b>Product Characteristics</b>	biologically inert physiologically harmless free of detectable phthalate acc. to REACH regulation
<b>Purity Grade</b>	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
<b>Sterilisation Process</b>	gas sterilisation by ethyleneoxid, 10-6
<b>Shelf - Life</b>	36 months at controlled storage conditions
<b>Packaging</b>	pcs per inner cs: 100 pcs per carton: 1.000 pcs per pallet: 70.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
<b>Approved by nerbe plus quality department</b>	