

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-812-8071	Revision: 2
Remarks		single use only

Specification

Product	np GREEN SWAB flocked swab, very fine nylon head, ABS stick, breakpoint at 80mm, sterile R, CE2797
Classification	CE2797
Dimensions	Length: 150 mm (+/- 3mm)
	Width: 2,5 mm (+/- 0,2mm)
	Diameter head: 3 mm (+/- 1mm)
Material Characteristics	head: nylon
	stick: Acrylonitrile butadiene styren (ABS)
Product Characteristics	biologically inert
	physiologically harmless
	free of detectable phthalate acc. to REACH regulation
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 100
	pcs per carton: 1.200
	pcs per pallet: 153.600
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

