

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-819-0110	Revision: 4
Remarks		single use only

Specification

Product	Wooden applicator with small cotton tip, Ø2,2x230mm, CE
Classification	CE
Dimensions	diameter: 2,20 mm (+/- 0,2 mm)
	length: 230,00 mm (+/- 0,2 mm)
	diameter tip: 5,00 - 5,50 mm (+/- 0,2 mm)
Material Characteristics	stick: wood (TI)
	tip: cotton (CO)
Product Characteristics	biologically inert
	physiologically harmless
	free of detectable phthalate acc. to REACH regulation
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Sterilisation Process	N/A
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 100
	pcs per carton: 10.000
	pcs per pallet: 360.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	