

Technical Datasheet

		EN ISO 13485 certified
Article - No:	09-819-8002	Revision: 5
Remarks		single use only

Specification

Product	Wooden applicator with small cotton tip, Ø2,2x150mm, single peel-pack, sterile EO, CE2797
Classification	CE2797
Dimensions	diameter stick: 2,20 mm (+/- 0,2 mm) length stick: 150,00 mm (+/- 0,2 mm)
Material Characteristics	stick: wood (TI) tip: cotton (CO)
Product Characteristics	biologically inert physiologically harmless free of detectable phthalate acc. to REACH regulation
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Sterilisation Process	gas sterilisation by ethyleneoxid, 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 100 pcs per carton: 1.000 pcs per pallet: 70.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	