

Technical Datasheet

		EN ISO 13485 certified
Article - No:	10-102-0100	Revision: 1
Remarks		single use only

Specification

Product	Micro test plate PP, 96 wells, 0,275ml, U-shape, 126,42x84,6x14,8 mm, transparent, np pcr ready, sterile R, max. RCF 4.000g, autoclavable
Dimensions	length: 126,42 mm (+/- 0,2 mm)
	width: 84,60 mm (+/- 0,2 mm)
	height: 14,80 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	0,275ml
Material Characteristics	polypropylene (PP)
Autoclavability	autocl.
Centrifugation	max. RCF 4.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 10
	pcs per carton: 100
	pcs per pallet: 6.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department

