

Technical Datasheet

		EN ISO 13485 certified
Article - No:	10-121-5100	Revision: 1
Remarks		single use only

Specification

Product	Micro test plate PS, 96 wells, 0,3ml, F-shape, SLAS/ANSI 1 & 4, highly transparent, np pcr ready, sterile R, max. RCF 4.500g
Dimensions	ANSI/SLAS 1 & 4
	length: 127,76 mm (+/- 0,5 mm)
	width: 85,48 mm (+/- 0,5 mm)
Colour	highly transparent
Nominal Capacity	0,3ml
Material Characteristics	polystyrene (PS)
Temperature range	-20°C - 60°C
Centrifugation	max. RCF 4.500g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 5
	pcs per carton: 100
	pcs per pallet: 6.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

